


Sostenil[®] *mini*

Sodium Hyaluronate



Restoring Synovial Balance in the Foot

Osteoarthritis of the foot

- Osteoarthritis (OA) of the foot is a painful condition that impacts on daily living, and is commonly associated with increasing age or with sports-related injuries.¹⁻³
- In general, OA of the foot is most common in the ankle, first tarsometatarsal and first metatarsophalangeal joints.⁴
- In particular, hallux rigidus or limitus, and hallux valgus severely reduce big toe flexibility causing pain when walking, running, dancing, etc.^{1,2}
- OA of the foot has received little clinical attention compared with OA of the knee or hip. More research is needed to evaluate the clinical impact of therapy on foot OA and to compare possible treatment options.
- Intra-articular injections of sodium hyaluronate are an effective treatment for OA of large joints, such as the knee and hip.⁵⁻⁸ In addition, they show promise for the treatment of OA in small joints, such as those in the foot.⁹⁻¹¹

Ostenil[®] mini: designed for relief from OA in small joints

- Ostenil[®] mini is indicated for the treatment of pain and restricted mobility associated with degenerative and traumatic changes in small synovial joints, such as those in the foot.
- Ostenil[®] mini is an isotonic solution of highly purified hyaluronan (sodium hyaluronate 1%) produced by bacterial fermentation. It is presented in individually packaged, sterile, pre-filled syringes, each containing 1 ml hyaluronan solution.
- Ostenil[®] mini is administered intra-articularly as a course of 1-3 injections at weekly intervals. The beneficial effects of treatment with Ostenil[®] mini can last up to 6 months.^{12,13}

A single injection of Ostenil[®] mini is an effective treatment for hallux limitus¹⁴

- The efficacy and tolerability of a single, intra-articular injection of Ostenil[®] mini were compared with those of steroid injections in a single-blind, randomized, 12-week study.
- A total of 37 patients with hallux limitus participated, and 40 joints were injected – 20 with Ostenil[®] mini and 20 with steroid.
- Within 2 weeks, Ostenil[®] mini significantly reduced pain on walking 20 m and the effects were maintained over 12 weeks ($p < 0.01$ vs baseline; Figure 1).




Figure 1. Ostenil[®] mini is significantly more effective than steroid at reducing pain on walking (TA=triamcinolone acetonide; VAS=visual analogue scale).

- Pain on walking was significantly reduced with Ostenil[®] mini at 4 and 8 weeks after injection compared with steroid ($p < 0.05$ vs steroid; Figure 1).
- Ostenil[®] mini and steroid significantly reduced pain on palpation or passive motion within 2 weeks of injection ($p < 0.01$ vs baseline).
- Compared with steroid, Ostenil[®] mini caused significant improvement in total score and pain evaluated on the Hallux Metatarsophalangeal-Interphalangeal Scale of the American Orthopaedic Foot and Ankle Society (AOFAS; $p < 0.05$; Figure 2).




Figure 2. Ostenil[®] mini is significantly more effective than steroid at improving total score on the AOFAS scale.

Multiple injections of Ostenil® mini are an effective treatment in hallux rigidus and hallux valgus¹⁵

- The efficacy and tolerability of a course of three intra-articular injections of Ostenil® mini were evaluated in a single-centre, exploratory, prospective, 15-week study.
- Ostenil® mini injections were administered at weekly intervals to 22 patients with Kellgren grades I-III hallux rigidus (n=11) or hallux valgus (n=11) and continuous pain for at least 3 months.
- Ostenil® mini significantly reduced pain, as measured on a visual analogue scale (VAS), within 1 week of the first injection ($p < 0.0001$ vs baseline). Pain reduction was maintained for at least 12 weeks after the last injection (Figure 3).
- Ostenil® mini reduced pain on walking and on load (Figure 4).
- Ostenil® mini significantly increased joint mobility ($p < 0.001$ vs pre-treatment).
- Overall, these changes resulted in improved walking ability (Figure 5).
- A total of 85.7% of patients reported recovery from initial OA symptoms at 15 weeks.
- Ostenil® mini was equally effective in treating hallux rigidus or hallux valgus.




Figure 3. Ostenil® mini significantly reduced pain within 1 week of injection, lasting for at least 12 weeks.




Figure 4. Ostenil® mini decreased pain on walking and on load.




Figure 5. Ostenil® mini improved walking ability.

Ostenil® mini is well tolerated

- Single injections of Ostenil® mini or steroid were equally well tolerated.¹⁴
- A course of three injections of Ostenil® mini was well tolerated, and no significant local or systemic adverse events were experienced.¹⁵

Ostenil[®] mini

Sodium Hyaluronate



- Ostenil[®] mini reduces pain associated with OA within 1 week of injection.¹⁵
- Ostenil[®] mini improves joint mobility and pain-free activity.^{14,15}
- A single injection of Ostenil[®] mini reduces pain and eases walking for at least 3 months.¹⁴
- A single injection of Ostenil[®] mini is at least as effective as an intra-articular injection of steroid.¹⁴
- Ostenil[®] mini is well tolerated.^{14,15}

References

1. Fajardo M and Di Cesare PE. *Drugs Aging* 2005; 22: 141–61.
2. Hockenbury RT. *Med Sci Sports Exerc* 1999; 31(7 Suppl.): S448–58.
3. Koh J and Dietz J. *Clin Sports Med* 2005; 24: 57–70.
4. McGuire JB. *Clin Podiatr Med Surg* 2003; 20: 469–85.
5. Dougados M et al. *Osteoarthritis Cartilage* 1993; 1: 97–113.
6. Altman RD and Moskowitz R. *J Rheumatol* 1998; 25: 2203–12.
7. Mathies B et al. Presented at the 5th Symposium of the International Cartilage Repair Society. May 26–29, 2004; poster 397.
8. Tikiz C et al. *Clin Rheumatol* 2005; 24: 244–50.
9. Marshall KW. *Foot Ankle Clin N Am* 2003; 8: 221–32.
10. Kelly M et al. *Am J Orthop* 2004; 33(Suppl. 2): 23–8.
11. Pleimann JH et al. *Foot Ankle Clin N Am* 2002; 7: 489–94.
12. Fuchs S et al. *Osteoarthritis Cartilage* 2006; 14: 82–8.
13. Bruhlmann P et al. Presented at the 10th World Congress on Osteoarthritis. December 8–11, 2005; poster P144.
14. Alvarez F et al. Presented at the 9th World Congress of the Osteoarthritis Research Society International. December 2–5, 2004; poster P347.
15. Talke M. *Arthritis + Rheuma* 2004; 24: 163–9 [article in German].

Prescribing information

OSTENIL[®] mini

Sodium hyaluronate 1.0%. Viscoelastic solution for injection into small joints.
Sterile by moist heat.

Composition: 1 ml isotonic solution contains 10.0 mg sodium hyaluronate from fermentation and sodium chloride, sodium monohydrogenphosphate, sodium dihydrogenphosphate, water for injection.

Indications: Pain and restricted mobility in degenerative and traumatic changes of small synovial joints, for example, the facet joints of the lumbar spine, the saddle joint of the thumb, the interphalangeal joints of the fingers and toes, the proximal joint of the big toe and the temporomandibular joint. In the treatment of bigger joints, for example the knee, hip or shoulder, OSTENIL[®] pre-filled syringes of 20 mg/2.0 ml should be used.

Contra-indications: OSTENIL[®] mini should not be used in patients with ascertained hypersensitivity to one of the constituents.

Precautions: Caution should be exercised in patients with known hypersensitivity to drugs. The general precautions for intra-articular injections should be observed, including measures to avoid joint infections. OSTENIL[®] mini should be injected accurately into the joint cavity, if necessary under imaging control. Avoid injections into blood vessels or surrounding tissues! As no clinical evidence is available on the use of hyaluronic acid in children, pregnant and lactating women or in inflammatory joint diseases such as rheumatoid arthritis or Bechterew disease, treatment with OSTENIL[®] mini is not recommended in these cases. Do not use if the pre-filled syringe or sterile pack are damaged. Store between 2°C and 25°C! Do not use after the expiry date indicated on the box. Keep out of the reach of children.

Side effects: Local secondary phenomena such as pain, feeling of heat, redness and swelling/joint effusion may occur in the joint treated with OSTENIL[®] mini. Application of an ice pack for five to ten minutes onto the treated joint will reduce the incidence of these events. In isolated cases the occurrence of side effects may be due to the intra-articular injection itself.

Interactions with other products: Avoid using OSTENIL[®] mini with instruments sterilised with quaternary ammonium salts solutions. No information on the incompatibility of OSTENIL[®] mini with other solutions for intra-

articular use is available to date. The concomitant use of an oral analgesic or anti-inflammatory drug during the first few days of treatment may be helpful for the patient.

Dosage and administration: Inject OSTENIL[®] mini into the affected joint once a week for a total of 1–3 injections. Several joints may be treated at the same time. Depending on the severity of the joint disease the beneficial effects of a treatment cycle may last at least six months. Repeat treatment cycles may be administered as required. In case of joint effusion it is advisable to reduce the effusion by aspiration, rest, application of an ice pack and/or intra-articular corticosteroid injection. Treatment with OSTENIL[®] mini can be started two to three days later.

The contents and outer surface of the OSTENIL[®] mini pre-filled syringe are sterile as long as the sterile pack remains unbroken. Take the pre-filled syringe out of the sterile pack, unscrew the Luer lock cap from the syringe, attach a suitable needle (for example 19 to 25 G) and secure it by turning slightly. If present remove the air bubble from the syringe prior to injection.

Characteristics and mode of action: Synovial fluid, which is viscoelastic due to the presence of hyaluronic acid, is found in all synovial joints, where it ensures normal, painless movement due to its lubricating and shock-absorbing properties. It is also responsible for the nutrition of the cartilage. In degenerative joint disorders such as osteoarthritis, the viscoelasticity of the synovial fluid is markedly reduced thereby decreasing its lubricating and shock-absorbing functions. This increases mechanical loading of the joint and cartilage destruction which ultimately results in pain and restricted mobility of the affected joint. Supplementing this synovial fluid with intra-articular injections of highly purified hyaluronic acid can ameliorate the viscoelastic properties of synovial fluid. This improves its lubricating and shock-absorbing functions and reduces mechanical overload of the joint. As a rule this results in a decrease in pain and an improvement in joint mobility which may last for several months after a treatment cycle.

Presentation and package size: One pre-filled syringe of 10 mg/1.0 ml OSTENIL[®] mini in a sterile pack.

OSTENIL[®] mini is a medical device. To be used by a physician only.

Last revision date: February 2006.

Manufacturer: TRB Chemedica AG, D-85540 Haar/München

TRB CHEMEDICA INTERNATIONAL SA
Geneva
Switzerland
Tel. + 41 22 703 49 00
Fax + 41 22 703 49 01
e-mail: international@trbchemedica.com
www.trbchemedica.com

Manufacturer:
TRB CHEMEDICA AG
Haar/Munich
Germany

